

AEON™



**INNOVATING
MINIMALLY INVASIVE
SURGERY WITH
S3 ENGINEERING™**



BOWA
MEDICAL

LEXINGTON
MEDICAL

AEON™

AEON WITH S₃ ENGINEERING™

Superior staple lines. Smooth articulation. Multi-speed gear.

Advanced reload architecture

Integrated I-beam assembly, fixed anvil and teardrop shape anvil buckets

High tech manufacturing

Proprietary manufacturing equipment incorporating AI

Cross compatibility

Cross compatible with Endo GIA™ and Signia™ powered handles



BETTER STAPLE FORMATION IN THICK TISSUE¹

AEON™ staple lines from a laparoscopic sleeve gastrectomy procedure



1. Data on file

SMOOTH RANGE OF MOTION WITH PRECISE ARTICULATION



Precise articulation

Freedom to articulate effortlessly for a more controlled firing

Multi-speed gear

Reduced trigger pull force for thicker tissue

Single-handed grasper

Quick, easy tissue grasping for increased efficiency



50%

LOWER FIRING FORCE WITH MULTI-SPEED GEAR¹ for controlled firing in thick tissue

1. Data on file

Reload Type	Cartridge Length	Product Code	Open Height	Closed Height	EA / BX
Grey	45mm	AESR45G	2.0mm	0.75mm	6
Grey	45mm	AESC45G	2.0mm	0.75mm	6
White	30mm	AESR30W	2.5mm	1.0mm	6
White Curved	30mm	ASR30WC	2.5mm	1.0mm	6
White	45mm	AESR45W	2.5mm	1.0mm	6
White Curved	45mm	AESC45W	2.5mm	1.0mm	6
Orange	45mm	AESR45R	3.25mm	1.5mm	6
Orange Curved	45mm	AESC45R	3.25mm	1.5mm	6
Purple	45mm	AESR45P	4.0mm	1.8mm	6
Purple Curved	45mm	AESC45P	4.0mm	1.8mm	6
White	60mm	AESR60W	2.5mm	0.75mm	6
White Short Tip	60mm	ASR60WS	2.5mm	1.0mm	6
White Curved	60mm	AESC60W	2.5mm	1.0mm	6
Orange	60mm	AESR60R	3.25mm	1.5mm	6
Orange Short Tip	60mm	ASR60RS	3.25mm	1.5mm	6
Orange Curved	60mm	AESC60R	3.25mm	1.5mm	6
Purple	60mm	AESR60P	4.0mm	1.8mm	6
Purple Curved	60mm	AESC60P	4.0mm	1.8mm	6
Purple Short Tip	60mm	ASR60PS	4.0mm	1.8mm	6
Black	60mm	AESR60B	5.0mm	2.2mm	6

Handle Type	Shaft Length	Product Code	EA / BX
Short	60mm	AESH060	3
Medium	160mm	AEHS160	3
Long	260mm	AEHS260	3

BOWA MEDICAL UK

 Bremridge House,
 Bremridge, Ashburton
 TQ13 7JX

T: 01364 652 426

F: 01364 653 589

E: support@bowa-medical.com

bowa-medical.co.uk