


# SHE SHA HANDLES



## The perfect mix – the SHE SHA smoke evacuation in combination with our length-adjustable smoke evacuation handles.


### NEW – SHE SHA handle with NON-Stick electrode

The handle extracts surgical smoke directly at their source. The new adjustable handle with a NON-STICK electrode provides excellent cutting and coagulation performance with even less sticking of tissue.


### NEW – More Flexibility

Both SHE SHA length-adjustable handles, with standard as well as with NON-Stick electrode, can be combined with our new needle electrode.

## SHE SHA HANDLES

		REF
	SHE SHA Handle, adjustable, knife, 3 m, single-use, sterile	802-033
	SHE SHA Handle, adjustable, NON-STICK knife, 3 m, single-use, sterile	802-034

## NEEDLE ELECTRODE (Shaft 2.4 mm, single-use, sterile)

		REF
	Needle electrode for SHE SHA handle 802-033 and 802-034	DNES620

For more information about our new smoke evacuation products or further SHE SHA products please visit our webshop: <https://int-shop.bowa-medical.com/en/Surgical-smoke/SHE-SHA/>

Contact us to request your personal demonstration. We are looking forward to hearing from you.

