




BOWA
M E D I C A L

1 **BOWA CABLE**
= 300
 **SINGLE-USE**
CABLES *

CABLES FOR ELECTROSURGERY

THE AWARD-WINNING COMBINATION OF DESIGN AND FUNCTIONALITY

BOWA MEDICAL cables are designed with the demands of the operating theatre in mind. With advanced safety features and a long service life of up to 300 sterilisation cycles, our reusable cables offer outstanding performance and value compared to single-use alternatives.



21p per use[†]

-  Monopolar and bipolar cables
-  Compatible with all common generators
-  Built to last
-  Ergonomic design
-  Advanced safety features



*Up to 300 sterilisation cycles per cable

[†]For our most popular bipolar cable purchased on NHS Supply Chain, if used 300 times

Award-winning design
thanks to ergonomically
perfected connectors

Optimum sealing
due to the cavity-free
2-component cast
design with smooth
surfaces

Safe connecting
orange marking shows
whether the plug is fully
connected to the socket

Safe identification
due to robust laser labelling

**Consistently safe
plug-in behaviour**
with stable spring
contacts

**Voltage safety
up to 6000 Vp**
thanks to double
insulation


Simple to use
thanks to the finger
depression trigger aid

Robust connectors
with kink protection
and strain relief

Highly flexible
due to the twist-free cable

**Clear identification and
protection against tripping**
with orange safety stripes

134 °C
Autoclave
273 °F



ONE-TIME OFFER

10% OFF

**PURCHASES OF 50+ CABLES
WITH THE CODE BOWA10**

Contact your local BOWA representative to find out more



LEARN MORE



FEEDBACK

BOWA
M E D I C A L